

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09257 - Goat Anti-TRPC7 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: TRPC7, transient receptor potential cation channel, subfamily C, member 7, TRP7, likely ortholog of mouse transient receptor potential cation channel, subfamily C, member 7, putative capacitative calcium channel

Official Symbol: TRPC7

Accession Number(s): NP_065122.1; NP_001161049.1; NP_001161048.1

Human GeneID(s): [57113](#)

Non-Human GeneID(s): 26946 (mouse), 282822 (rat)

Important Comments: This antibody is expected to cross-react with isoform 1, 2 and 3 (NP_065122.1; NP_001161049.1; NP_001161048.1).

Immunogen

Peptide with sequence KLCKSKAKSCEND, from the internal region of the protein sequence according to NP_065122.1; NP_001161049.1; NP_001161048.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: A 95kDa band has been observed in Chicken Embryo Heart lysates in the following paper: Sabourin et al, Cardiovasc Res. 2011 Jun 13, PMID 21672930.

Immunoprecipitation: This antibody has been successfully used in IP on lysates of Chicken Embryo Heart (PMID: 21672930).

Species Reactivity

Tested: Chicken

Expected from sequence similarity: Human, Mouse, Rat, Dog

Specific References

This antibody has been successfully used (in WB on Chicken heart (atrial and ventricular myocytes)) in the following paper:

Sabourin J, Antigny F, Robin E, Frieden M, Raddatz E.

Activation of Transient Receptor Potential Canonical 3 (TRPC3)-mediated Ca²⁺ Entry by A1 Adenosine Receptor in Cardiomyocytes Disturbs Atrioventricular Conduction.

J Biol Chem. 2012 Aug 3;287(32):26688-701.

PMID: 22692208

This antibody has been successfully used in the following paper:

Sabourin J, Robin E, Raddatz E.

A key role of TRPC channels in the regulation of electromechanical activity of the developing heart.

Cardiovasc Res. 2011 Jun 13. [Epub ahead of print]

PMID: 21672930