

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB09658 - Goat Anti-UBR2 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: UBR2, ubiquitin protein ligase E3 component n-recognin 2, C6orf133, DKFZp686C08114, KIAA0349, MGC71112, RP3-392M17.3, bA49A4.1, dJ242G1.1, dJ392M17.3, OTTHUMP00000016403, OTTHUMP00000016405, OTTHUMP00000039767, OTTHUMP00000039768, ubiquitin ligase E3 alpha-II

Official Symbol: UBR2

Accession Number(s): NP_056070.1

Human GeneID(s): [23304](#)

Non-Human GeneID(s): 224826 (mouse)

Immunogen

Peptide with sequence SCPKSGGDKSRAPT, from the internal region of the protein sequence according to NP_056070.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:32000.

Western blot: Preliminary experiments gave bands at approx 150kDa in Rat Testis lysates after 0.01µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the bands we observe given the calculated size of 199kDa according to Rat NP_001171542.1. The detected band was successfully blocked by incubation with the immunizing peptide (and BLAST results with the immunizing peptide sequence did not identify any other proteins to explain the additional bands).

IHC: Paraffin embedded Brain (Cortex). Recommended concentration: 5µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Dog, Cow

EB09658 (5µg/ml) staining of paraffin embedded Human Cortex. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.