

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB06994 - Goat Anti-UCP1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: UCP1, uncoupling protein 1 (mitochondrial, proton carrier), HGNC:12517, SLC25A7, UCP, mitochondrial brown fat uncoupling protein, uncoupling protein 1

Official Symbol: UCP1

Accession Number(s): NP_068605.1

Human GeneID(s): [7350](#)

Immunogen

Peptide with sequence C-EQLKRELSKSRQ, from the C Terminus of the protein sequence according to NP_068605.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:128000.

Western blot: Approx 35kDa band observed in Human Adipose lysates (calculated MW of 33.0kDa according to NP_068605.1). Recommended concentration: 0.1-1 µg/ml. Primary incubation was 1 hour.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog

EB06994 (1µg/ml) staining of Human Adipose lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.