

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05991 - Goat Anti-USP15 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: USP15, ubiquitin specific protease 15, UNPH4, KIAA0529, deubiquitinating enzyme, ubiquitin specific peptidase 15, MGC131982, MGC149838, MGC74854

Official Symbol: USP15

Accession Number(s): NP_006304.1

Human GeneID(s): [9958](#)

Immunogen

Peptide with sequence NDIENENCMHTN, from the C Terminus of the protein sequence according to NP_006304.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:2000.

IHC: In paraffin embedded Human Liver shows patterned intracellular staining in hepatocytes. Recommended concentration, 2.5-3.8µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Rat, Dog, Pig, Cow



EB05991 (2µg/ml) staining of paraffin embedded Human Liver. Steamed antigen retrieval with citrate buffer pH 6, HRP-staining.