

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05972 - Goat Anti-USP4 / UNP Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: USP4, ubiquitin specific protease 4 (proto-oncogene), UNP, Unph, ubiquitin specific protease, proto-oncogene, ubiquitin specific peptidase 4 (proto-oncogene), MGC149848, MGC149849, deubiquitinating enzyme 4, ubiquitin carboxyl-terminal esterase 4, ubiquitin carboxyl-terminal hydrolase 4, ubiquitin specific protease 4, ubiquitin thiolesterase 4, ubiquitin-specific processing protease 4

Official Symbol: USP4

Accession Number(s): NP_003354.2; NP_955475.1

Human GeneID(s): [7375](#)

Immunogen

Peptide with sequence AEGGGCRERPDAE, from the N Terminus of the protein sequence according to NP_003354.2; NP_955475.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Western blot: Approx 100kDa band observed in Human Heart lysates (calculated MW of 104kDa according to NP_955475.1). Recommended concentration: 1-3µg/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



Approx 100kDa band observed in Human Heart lysates (calculated MW of 104kDa according to NP_955475.1).
Recommended concentration: 1-3µg/ml.