

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre
77 Heyford Park
Upper Heyford
Oxfordshire
OX25 5HD
UK

Enquiries:

info@everestbiotech.com

Sales:

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

**Research Use Only. Not for
diagnostic or therapeutic use.**

EB05974 - Goat Anti-USP6 / TRE2 / TRE17 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: USP6, ubiquitin specific protease 6 (Tre-2 oncogene), HRP1, TRE2, Tre2, TRE17, Tre-2, tre-2 oncogene, hyperpolymorphic gene 1, USP6-short, ubiquitin specific peptidase 6-, ubiquitin specific protease 6, ubiquitin specific peptidase 6 (Tre-2 oncogene), ubiquitin-specific protease USP6

Official Symbol: USP6

Accession Number(s): NP_004496.1

Human GeneID(s): [9098](#)

Immunogen

Peptide with sequence KISPLHHLQMECSP, from the C Terminus of the protein sequence according to NP_004496.1.

Please note the [peptide](#) is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested


Peptide ELISA: antibody detection limit dilution 1:16000.

Western blot: Approx 90kDa band observed in Human Placenta lysates (calculated MW of 89.5kDa according to NP_004496.1). Recommended concentration: 1-3µg/ml. An additional band of 50kDa was consistently observed, however this band was not blocked by the immunizing peptide and it is therefore a non-specific signal. We call for caution when used for other assays than Western blot.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human



EB05974 (2 μ g/ml) staining of Human Placenta lysate (35 μ g protein in RIPA buffer) with (B) and without (A) blocking with the immunizing peptide. Primary incubation was 1 hour. Detected by chemiluminescence.